

*Newheek*

---

# X-Ray Hand Switch

## -----Type L01

*Weifang Newheek Medical Image Equipment Co., Ltd.*

Add: No. 1 Factory Building of No.13426 of Yuqing East Street,  
Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



**Weifang Newheek Medical Image Equipment Co., Ltd.**

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

## Quick Details

Properties: Medical X-ray Equipments & Accessories

Brand Name: Newheek

Product Name: X-ray exposure hand switch

Model Number: L01

Place of Origin: Shandong, China (Mainland)

Application: x ray machine

2 steps design : step 1 exposure ready status, step 2 control exposure

Current: 3A

Housing materials: ABS plastic

Operation voltage: 125v

Standard core: 3 core or 4 core

Standard Length: 3m wire (30cm spring coil wire)

Mechanical life: 1,000,000 times

Electric life: 200,000 times

Structure: omron inner switch

Customization: Available

## Product Application

### **Applicable to x-ray exposure of radiography or fluoroscopy equipment**

X ray hand exposure hand switch is mainly used on portable x ray, mobile x ray, stationary x ray, analog x ray, digital x ray, radiography x ray etc x ray equipment.

As well as it is applicable to beauty laser device, healthy recovery device etc field.

## HighLight

- Streamline design outlook, comfortable handle.
- Housing adopt ABC engineering plastic, anti-shock and fire retardant.
- Inner adopt OMRON micro switch, high sensitivity.
- Mechanical life of 1,000,000 times above.
- Electrical life of 200,000 times above.

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

Detail diagram



Picture of Hand switch L01



Picture of Handle and cradle

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



Picture of **Two-Button**



Picture of **3m-wire** (30cm spring coil wire)

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

### Technical Features

#### 3-core switch

Operation Voltage	Current	Housing Material	Hand Switch Cable Core		
			Red	White	Green
220V AC	10A	white, ABS plastic	I Step	Common Point	II Step
Ambient Temperature		Relative Humidity		Atmospheric Pressure	
-20°C-70°C		≤ 93%		50-106kPa	

#### 4-core switch

Operation Voltage	Current	Housing Material	Hand Switch Cable Core			
			White	Black	Red	Green
220V AC	10A	white, ABS plastic	I Step		II Step	
Ambient Temperature		Relative Humidity		Atmospheric Pressure		
-20°C-70°C		≤ 93%		50-106kPa		

### Configuration:

Standard	3-core, 3m-wire	4-core, 3m-wire
Customized	3-core, 6m or 10m wire	4-core, 6m or 10m wire
	Other requirement	Other requirement

**Weifang Newheek Medical Image Equipment Co., Ltd.**

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)

## Our Service

**Original manufacturer of x-ray machine accessories hand switch and foot switch for more than 16 years.**

- √ Customers could find all kinds of x-ray machine parts here.
- √ Offer on line technological support.
- √ Promise super product quality with best price and service.
- √ Support the third part inspection before delivery.
- √ Ensure the shortest delivery time.

## Contact us

Email : [info@xrayswitch.com](mailto:info@xrayswitch.com)  
Hot Line: +86 536 8882360

## FAQ

The common problems on the **hand switch L01**

1. Can you accept customize hand switch?

Yes. Please feel free to send us your detail requirement by email here.  
You will get our reply in 24 hours.

2. Can we print our logo/website/company name on products?

Yes, please advise the size and Pantone Code of the logo.

3. What's the lead time for regular order?

3-5 days for order quantity lower than 100pc hand switch; For more quantity, usually 15 days.

4. Can I get a discount?

Yes, for the order quantity more than 50 pcs hand switch, please contact us to get the best price.

5. What is the quality warranty?

We offer 12 months warranty after the hand switch is delivered.

***Weifang Newheek Medical Image Equipment Co., Ltd.***

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: [www.xrayswitch.com](http://www.xrayswitch.com)

Tel: 86 536 8882360

Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



Certificate EMC

Certificate of Conformity

Certificate Number: B-E16029006  
EMC Directive 2014/30/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.  
Address.....: The North 2nd Floor of No.1 Factory Building of No.13426 Yuqing East Street, Hi-Tech Industrial Development Zone, Weifang, Shandong, China  
Manufacturer.....: Same As Holder  
Product.....: Hand Switch  
Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06, L07, L08, L09, L10  
Technical Data.....: AC 110-240V~, 10A / DC 30V, 5A

The submitted products have been tested by us with the following standard(s) and found to be in compliance with the listed European Directives.

- EN 61000-6-3:2007+A1:2011;
- EN 61000-6-1:2007;
- EN 61000-3-2:2006+A1:2009+A2:2009;
- EN 61000-3-3:2013

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.



*[Signature]*

Certification Manager  
Date: Feb. 24, 2016



Company No.07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indecon Square, London, United Kingdom  
China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China  
Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.  
Website: [www.xrayswitch.com](http://www.xrayswitch.com) Tel: 86 536 8882360 Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)



Certificate ROHS

Certificate of Conformity

Certificate Number: B-R16029007  
RoHS Directive 2011/65/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.

Address.....: The North 2nd Floor of No.1 Factory Building of No.13426  
Yuqing East Street, Hi-Tech Industrial Development Zone,  
Weifang, Shandong, China

Manufacturer.....: Same As Holder

Product.....: Hand Switch

Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06,  
L07, L08, L09, L10

The submitted products have been tested by us and found to be in compliance with the listed European Directives.

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.



Certification Manager  
Date: Feb. 24, 2016



Company No. 07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indecon Square, London, United Kingdom  
China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China  
Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.  
Website: [www.xrayswitch.com](http://www.xrayswitch.com) Tel: 86 536 8882360 Email: [info@xrayswitch.com](mailto:info@xrayswitch.com)