

Nöwhëëk

X Ray Exposure Hand Switch

—— Type L04

Website : www.xrayswitch.com

Email : info@xrayswitch.com

Service Line : +86 536 8882360

Quick Details

Product Name: X-ray exposure hand switch
Place of Origin: Shandong, China (Mainland)
Brand Name: Newheek
Model Number: NKL04
Protection Level: IP67
Max. Current: 10A
Max. Voltage: 220VAC
Product name: two step hand switch
Inner switch: OMRON microswitch
Design: double knobs double microswitch
Type: 2-steps
Electrical Life: 300,000 times
Mechanical Life: 50,000,000 times
Housing material: ABC plastic
Length: 3m, 6m, 10m
Standard: 3 core, 4 core
Customization: Available

Product Application

Applicable to x-ray exposure of radiography or fluoroscopy equipment

X ray hand exposure hand switch is mainly used on portable x ray, mobile x ray, stationary x ray, analog x ray, digital x ray, radiography x ray etc x ray equipment.

As well as it is applicable to beauty laser device, healthy recovery device etc field.

HighLight

- Streamline design outlook, comfortable handle.
- Housing adopt ABC engineering plastic, anti-shock and fire retardant.
- Inner adopt OMRON micro switch, high sensitivity.
- Mechanical life of 50,000,000 times above.
- Electrical life of 300,000 times above.

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Detail diagram



Picture of Hand switch L04



Picture of Handle and cradle

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com



Picture of **Handle**



Picture of **3m-wire** (30cm spring coil wire)

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Technical Features

3-core switch

Operation Voltage	Current	Housing Material	Hand Switch Cable Core		
			Red	White	Green
220V AC	10A	white, ABS plastic	I Step	Common Point	II Step
Ambient Temperature		Relative Humidity		Atmospheric Pressure	
-20°C-70°C		≤ 93%		50-106kPa	

4-core switch

Operation Voltage	Current	Housing Material	Hand Switch Cable Core			
			White	Black	Red	Green
220V AC	10A	white, ABS plastic	I Step		II Step	
Ambient Temperature		Relative Humidity		Atmospheric Pressure		
-20°C-70°C		≤ 93%		50-106kPa		

Configuration:

Standard	3-core, 3m-wire	4-core, 3m-wire
Customized	3-core, 6m or 10m wire	4-core, 6m or 10m wire
	Other requirement	Other requirement

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Our Service

Original manufacturer of x-ray machine accessories hand switch and foot switch for more than 16 years.

- √ Customers could find all kinds of x-ray machine parts here.
- √ Offer on line technological support.
- √ Promise super product quality with best price and service.
- √ Support the third part inspection before delivery.
- √ Ensure the shortest delivery time.

Contact us

Email : info@xrayswitch.com

Hot Line: +86 536 8882360

FAQ

The common problems on the **hand switch L04**

1. Can you accept customize hand switch?

Yes. Please feel free to send us your detail requirement by email here.
You will get our reply in 24 hours.

2. Can we print our logo/website/company name on products?

Yes, please advise the size and Pantone Code of the logo.

3. What's the lead time for regular order?

3-5 days for order quantity lower than 100pc hand switch; For more quantity, usually 15 days.

4. Can I get a discount?

Yes, for the order quantity more than 50 pcs hand switch, please contact us to get the best price.

5. What is the quality warranty?

We offer 12 months warranty after the hand switch is delivered.

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Certificate EMC

Certificate of Conformity

Certificate Number: B-E16029006

EMC Directive 2014/30/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.

Address.....: The North 2nd Floor of No.1 Factory Building of No.13426 Yuqing East Street, Hi-Tech Industrial Development Zone, Weifang, Shandong, China

Manufacturer.....: Same As Holder

Product.....: Hand Switch

Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06, L07, L08, L09, L10

Technical Data.....: AC 110-240V~, 10A / DC 30V, 5A

The submitted products have been tested by us with the following standard(s) and found to be in compliance with the listed European Directives.

- EN 61000-6-3:2007+A1:2011;
- EN 61000-6-1:2007;
- EN 61000-3-2:2006+A1:2009+A2:2009;
- EN 61000-3-3:2013

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.



[Signature]

Certification Manager
Date: Feb. 24, 2016



Company No.07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indecon Square, London, United Kingdom
 China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China
 Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
 Website: www.xrayswitch.com Tel: 86 536 8882360 Email: info@xrayswitch.com

Certificate ROHS

Certificate of Conformity

Certificate Number: B-R16029007
RoHS Directive 2011/65/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.

Address.....: The North 2nd Floor of No.1 Factory Building of No.13426 Yuqing East Street, Hi-Tech Industrial Development Zone, Weifang, Shandong, China

Manufacturer.....: Same As Holder

Product.....: Hand Switch


Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06, L07, L08, L09, L10

The submitted products have been tested by us and found to be in compliance with the listed European Directives.

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.




 Certification Manager
 Date: Feb. 24, 2016



Company No. 07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indecon Square, London, United Kingdom
 China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China
 Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
 Website: www.xrayswitch.com Tel: 86 536 8882360 Email: info@xrayswitch.com