

Newhëëk

X Ray Hand Switch

——— Type L06

Website: www.xrayswitch.com

Email : info@xrayswitch.com

Service line : +86 536 8882360

Quick Details

Properties: Medical X-ray Equipments & Accessories
Brand Name: Newheek
Product Name: X-ray exposure hand switch
Model Number: L06
Place of Origin: Shandong, China (Mainland)
Application: x ray machine
Type : step 1 control exposure
Function :One step to control dental x ray expose)
LED : No
Current: 3A
Housing materials: ABS plastic
Operation voltage: 125v
Standard core: 2 core
Length: 3m 6m 10m
Mechanical life: 1,000,000 times
Electric life: 200,000 times
Structure: omron inner switch
Customization:Available

Product Application

Applicable to x-ray exposure of radiography or fluoroscopy equipment

X ray hand exposure hand switch is mainly used on portable x ray, mobile x ray, stationary x ray, analog x ray, digital x ray, radiography x ray etc x ray equipment.

As well as it is applicable to beauty laser device, healthy recovery device etc field.

HighLight

- Streamline design outlook, comfortable handle.
- Housing adopt ABC engineering plastic, anti-shock and fire retardant.
- Inner adopt OMRON micro switch, high sensitivity.
- Mechanical life of 1,000,000 times above.
- Electrical life of 200,000 times above.

Weifang Newheek Medical Image Equipment Co., Ltd.


Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Detail diagram



Picture of **Type L06 X-ray Hand Switch**



Picture of **One-step Button (No LED)**


Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.


Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com



Picture of **2-core Wire**



Picture of **6m-Wire** (60cm spring coil wire)

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Technical Features

Operation Voltage	Current	Housing Material	Hand Switch Cable Core	
			Red	Green
125V	3A	white, ABS plastic	I Step	
Ambient Temperature		Relative Humidity	Atmospheric Pressure	LED
-20°C-70°C		≤ 93%	50-106kPa	no

Technical Features

Standard	2-core, 3m-wire
Customized	2-core, 6m or 10m wire
	Other requirement

Our Service

Original manufacturer of x- ray machine accessories hand switch and foot switch for more than 16 years.

- √ Customers could find all kinds of x-ray machine parts here.
- √ Offer on line technological support.
- √ Promise super product quality with best price and service.
- √ Support the third part inspection before delivery.
- √ Ensure the shortest delivery time.

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Contact us

Email : info@xrayswitch.com

Hot Line: +86 536 8882360

The common problems on the **hand switch L06**

FAQ

1. Can you accept customize hand switch?

Yes. Please feel free to send us your detail requirement by email here.
You will get our reply in 24 hours.

2. Can we print our logo/website/company name on products?

Yes, please advise the size and Pantone Code of the logo.

3. What's the lead time for regular order?

3-5 days for order quantity lower than 100pc hand switch; For more quantity, usually 15 days.

4. Can I get a discount?

Yes, for the order quantity more than 50 pcs hand switch, please contact us to get the best price.

5. What is the quality warranty?

We offer 12 months warranty after the hand switch is delivered.

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.

Website: www.xrayswitch.com

Tel: 86 536 8882360

Email: info@xrayswitch.com

Certificate EMC

Certificate of Conformity

Certificate Number: B-E16029006
EMC Directive 2014/30/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.
Address.....: The North 2nd Floor of No.1 Factory Building of No.13426 Yuqing East Street, Hi-Tech Industrial Development Zone, Weifang, Shandong, China
Manufacturer.....: Same As Holder
Product.....: Hand Switch
Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06, L07, L08, L09, L10
Technical Data.....: AC 110-240V~, 10A / DC 30V, 5A

The submitted products have been tested by us with the following standard(s) and found to be in compliance with the listed European Directives.

- EN 61000-6-3:2007+A1:2011;
- EN 61000-6-1:2007;
- EN 61000-3-2:2006+A1:2009+A2:2009;
- EN 61000-3-3:2013

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.



[Signature]

Certification Manager
Date: Feb. 24, 2016



Company No.07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indecon Square, London, United Kingdom
China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China
Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xrayswitch.com Tel: 86 536 8882360 Email: info@xrayswitch.com

Certificate ROHS

certificate of conformity

Certificate Number: B-R16029007
RoHS Directive 2011/65/EU



Holder.....: Weifang Newheek Electronic Technology Co.,Ltd.

Address.....: The North 2nd Floor of No.1 Factory Building of No.13426 Yuqing East Street, Hi-Tech Industrial Development Zone, Weifang, Shandong, China

Manufacturer.....: Same As Holder

Product.....: Hand Switch


Model No.....: XC-01, XC-02, XC-03, XC-04, L01, L02, L03, L04, L05, L06, L07, L08, L09, L10

The submitted products have been tested by us and found to be in compliance with the listed European Directives.

The test results apply only to the particular sample tested and to the specific tests carried out. Technical Report and documentation are at the Holder's disposal.

This certificate applies specifically to the sample investigated in our test reference number only. The CE markings as shown below can be affixed on the product after preparation of necessary technical documentation. Other relevant Directives have to be observed.




 Certification Manager
 Date: Feb. 24, 2016



Company No. 07113834



Beide (UK) Product Service Limited

U.K.: Flat 107, 25 Indecon Square, London, United Kingdom
 China: 6F, Bldg E, Hourui 3rd Ind Zone, Xixiang, Bao'an Dist, Shenzhen, China
 Http://www.szbeide.com E-mail: admin@szbeide.com

Weifang Newheek Medical Image Equipment Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
 Website: www.xrayswitch.com Tel: 86 536 8882360 Email: info@xrayswitch.com